

**Safety Communication**

**رسالة سلامة**

**Use of an Isolation Transformer in combination with Sentec's Digital Monitor  
in home care settings**

<b>Device/ Product Description:</b>	Digital Monitor	
<b>Ref. Numbers:</b>	SDM	
<b>Manufacturer:</b>	SenTec AG.	
<b>Problem:</b>	Non-compliance of the affected devices with normative requirements from IEC 60601-1-11:2015 Medical electrical equipment when used in home environments.	
<b>Recommendation /Actions:</b>	<ul style="list-style-type: none"> <li>- Make sure this document is reached to the end-user</li> <li>- Use of an isolation transformer in combination with the SDM in home care settings to comply with applicable safety standards.</li> </ul> <p>For more information, please click <a href="#">here</a>. If you think you had a problem with your device or a device your patient uses, please report the problem to SFDA through: <a href="#">NCMDR</a> <a href="#">Vigilance system</a> (19999) unified call center</p>	
<b>Authorized Representative Details</b>	AR name:	Medical Regulations Gate
	Assigned Contact Person:	Abdullah Alsaiani
	Mobile/Phone:	+966 540640847
	Email:	ameen@mr-gate.com